

Dermalux Indication Protocol Matrix

SKIN CONCERN	KEY INDICATIONS	DERMALUX PROTOCOL	PROGRAMME	WAVELENGTHS		
Skin Ageing	Fine lines and wrinkles, loss of elasticity, enlarged pores, puffiness	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with anti-ageing treatments	Skin Ageing Skin Ageing Post Treatment Rejuvenation		●	●
Skin Ageing: Pigmentation	Sun damage, age spots, uneven skin tone	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with brightening treatments	Pigmentation Pigmentation Post Treatment Rejuvenation		●	●
Skin Ageing: Dry skin	Dry and tight, rough skin texture, flaking skin	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with resurfacing and soothing treatments	Skin Ageing Skin Ageing Post Treatment Rejuvenation		●	●
Skin Ageing: Complexion	Dull sallow skin, dehydration, lack of radiance	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with resurfacing and hydrating treatments	Skin Ageing Skin Ageing Post Treatment Rejuvenation		●	●
Pigmentation	Hormonal pigmentation, post inflammatory pigmentation, photoallergic and phototoxic responses	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with brightening treatments	Pigmentation Pigmentation Post Treatment Rejuvenation		●	●
Spot Prone Skin: Mild to Moderate	Occasional blemishes, comedones, combination skin	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with mild spot treatments	Spot Prone Skin Spot Prone Skin Post Treatment Anti-Bacterial	●	●	●
Spot Prone Skin: Moderate to Severe	Persistent blemishes and breakouts, enlarged and sore spots, oily skin	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with advanced spot treatments	Spot Prone Skin Spot Prone Skin Post Treatment Anti-Bacterial	●	●	●



SKIN CONCERN	KEY INDICATIONS	DERMALUX PROTOCOL	PROGRAMME	WAVELENGTHS		
Redness: Vascular	Visible capillaries, thin and fragile skin	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with vascular and calming treatments	Redness Redness Post Treatment Calming		● ●	● ● ●
Redness: Skin Tone	Red skin tone, flushed complexion, diffused redness	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with calming treatments	Redness Redness Post Treatment Calming		● ●	● ● ●
Sensitive Skin: Irritation	Compromised, irritated, reactive skin	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments ENHANCE: Combine with calming treatments	Sensitive Skin Sensitive Skin Post Treatment Calming			● ● ●
Sensitive Skin: Problem conditions	General sensitive and fragile skin conditions	CORRECT: Dermalux Power Plan MAINTAIN: Dermalux Booster treatments	Tri-Wave Therapy Tri-Wave Therapy	● ●	● ●	● ●

DERMALUX POWER PLAN	8 week intensive programme. 12 treatments; 2 treatments per week for 4 weeks followed by 1 treatment per week for 4 weeks
DERMALUX BOOSTER TREATMENTS	Individual treatments. Recommend monthly or as required
DERMALUX ENHANCER TREATMENTS	Apply immediately post treatment for superior results